



**SAN MATEO COUNTY PUBLIC HEALTH LABORATORY
TEST INFORMATION**

Section: Syphilis Serology

Test Name: *Treponema pallidum* – Particle Agglutination (TP-PA)

Test Includes: <i>Treponema pallidum</i> – Particle Agglutination	
<i>Reporting</i>	
Results Available: 7 days from receipt	Contact Number: (650) 573-2500
<i>Reference</i>	
Method: Agglutination	
Turnaround Time: 7 days	Reference Range: Non-Reactive
Limitations: May not detect a recent infection, or infection in a person with a severely compromised immune system.	Interpretation: Non-reactive indicates that the patient does not have detectable antibody to the infectious agent. Reactive indicates that the patient has detectable antibody to the infectious agent, and depending on the clinical picture, may have a current or past infection.
<i>Specimen Requirements</i>	
Specimen Collection: Venipuncture	Sample Type: Blood
Volume/Amount Required: 7-10 mls whole blood or 5 ml serum or plasma	Preferred Specimen: Serum, plasma
Collection/Preservation: Red top or tiger top tube	Storage Instructions: Do not freeze or refrigerate whole blood. Separated serum may be held at 2-8° C.
Causes for Rejection: Discrepancy in specimen identification; insufficient quantity of specimen; gross hemolysis	Sample Container: Red top or tiger top tube
Sample Test Kit:	Availability: Performed Friday
<p>Diagnostic Information: The Serodia TP-PA is a qualitative gelatin particle agglutination assay for the detection of <i>Treponema pallidum</i> antibodies in human serum. The test is based on the agglutination of colored gelatin particles sensitized with <i>T. pallidum</i> antigen. Serum containing specific antibodies will react with the antigen sensitized gelatin particles to form a smooth mat of agglutinated particles. A compact button formed by the settling of non-agglutinated particles characterizes negative reactions. Unlike the nontreponemal tests, once the TP-PA test becomes reactive, it remains reactive for many years. Since the reactivity found with the TP-PA does not indicate response to therapy, it is not suitable for monitoring treatment. The TP-PA test does not distinguish between syphilis and other treponematoses such as yaws, pinta and bjel. TP-PA is run only if the RPR card test results are reactive.</p>	
<i>Specimen Submission</i>	
Request Form: Standard Clinical Test Request form or through computer interface	Specimen Handling: Use Universal Precautions
Transport Temperature: Ambient temperature for clotted blood; separated serum	Shipping Requirements: Ship on cold packs (2-8°C), or frozen (-20°C) and mailed on dry



**SAN MATEO COUNTY PUBLIC HEALTH LABORATORY
TEST INFORMATION**

at 2-8°C (refrigerated) or -20°C (frozen).	ice.
<i>Billing</i>	
CPT Code(s): 86781	Fees: \$15.00
Effective Date: July 1, 2008	