


Live 06/17/2015

This is a new update in Soarian wherein MD or other providers can now view the digital radiology images on the Patient Record without logging on to the Picom Web.

There are 2 ways to view patient's radiology result images;

1. From the Patient Record – click on the icon  and the result image will open.



TEST, ADAM | DOB: 02/09/2004(11y) ♀ | MRN: 1005255 | ACCT#: 20218545 | 39TH AVE CLINIC | Admit Date: 05/26/2015


All Occurrences | 02/26/2015 to 05/27/2015 | Last 24 Hours

Category	Procedure	Date/Time	Image Icon	Report Status
Dx Radiology	Foot Ap & Lat Lt	05/11/15 15:23		400135***F ina
	Hand And Wrist Lt	12/08/14 10:46		***Final Repo
	Low Extr Entire Child 2Vw Lt	06/24/14 16:09		***Final Repo...
CT Scan	CT Abdomen W Con	05/11/15 15:39		400125***F ina
MRI	MRI-Chest-W Con	03/06/15 16:46		***Final Repo
	MRI-Head-W Con	10/09/14 09:51		***Final Repo

SMMC_SOARIAN	TITLE: HOW TO VIEW RADIOLOGY IMAGE IN SOARIAN		
EXTERNAL	VERSION: 1.0	DATE: 06/15/2015	PREPARED: CECIL AGDIPA
PAGE 1 OF 3	DOCUMENT REF: SPT/TP-06-15, V.1.0		AUTHORIZED: CECIL AGDIPA

How to View Radiology Image in Soarian

2. From the Patient Text Result , click on the View image button, the result will open.




Click on the image series on the left and the image will show on the right side of the screen



SMMC_SOARIAN	TITLE: HOW TO VIEW RADIOLOGY IMAGE IN SOARIAN		
EXTERNAL	VERSION: 1.0	DATE: 06/15/2015	PREPARED: CECIL AGDIPA
PAGE 2 OF 3	DOCUMENT REF: SPT/TP-06-15, V.1.0		AUTHORIZED: CECIL AGDIPA

How to View Radiology Image in Soarian

Scroll the middle with your mouse to view the different angle of the images.



SMMC_SOARIAN	TITLE: HOW TO VIEW RADIOLOGY IMAGE IN SOARIAN		
EXTERNAL	VERSION: 1.0	DATE: 06/15/2015	PREPARED: CECIL AGDIPA
PAGE 3 OF 3	DOCUMENT REF: SPT/TP-06-15, V.1.0		AUTHORIZED: CECIL AGDIPA