

Section: Environmental Test Name: Heterotrophic Plate Count

Test Includes: Detection of heterotrophic bacteria	
Reporting	
Results Available: Wednesday through Friday	Contact Number : 650-573-2500
Reference	
Method: Standard plate count	
Turnaround Time: 2-3 days	Reference Range : Varies by sample type
Limitations : No Thursday or Friday deliveries. Specimen must delivered within 4 hours of	Interpretation: Varies by sample type
collection.	
Specimen Requirements	
Specimen Collection:	Sample Type : Pool, sauna, hot tub or potable
	water
Volume/Amount Required: 25-50 ml	Preferred Specimen:
Collection/Preservation: sodium thiosulfate	Storage Instructions : 2-4° C
Causes for Rejection: Too long in transit;	Sample Container: Approved 120 ml bottle
sample not refrigerated	
Sample Test Kit:	Availability: Monday through Wednesday
	only. No emergency samples from the public.
Test Information : Detection limits – 0-3,000 colonies/ml; extra dilutions on request.	
Specimen Submission	
Request Form: Environmental Test Request	Specimen Handling:
Transport Temperature : <10 ° C with freezer	Shipping Requirements:
packs	
Billing	
Effective Date: June 15, 2011	Fees : \$25.00